



Recombinant Rat Proteinase-activated receptor 1 (F2r)

Product Code	CSB-CF007924RA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P26824
Product Type	Transmembrane Protein
Immunogen Species	Rattus norvegicus (Rat)
Sequence	SFFLRNPSEDTFEQFPLGDEEEKNESIPLEGRAVYLNKSRFPPMPPPPFISEDA SGYLTS PWLTLFIPSVYTFVFIVSLPLNILAIAVFVFRMKVKKPAVVYMLHLAMADVLFVSV LPFK ISYYFSGTDWQFGSGMCRFATAACYCNMYASIMLMTVISIDRFLAVVYPIQSL WRTLGR ANFTCVVIWVMAIMGVVPLLLKEQTTQVPGLNITTCHDVLNETLLHGFYSYYFS AFSAIF FLVPLIISTVCYTSIIRCLSSSAVANRSKKSRAFLSAAVFCIFIVCFGPTNVLLIVH YL LLSDSPGTETAYFAYLLCVCVTSVASCIDPLIYYYASSECCQKHLYSILCCRESSD SNSCN STGQLMPSKMDTCCSHLNNSIYKLLA
Source	in vitro E.coli expression system
Target Names	F2r
Protein Names	Recommended name: Proteinase-activated receptor 1 Short name= PAR-1 Alternative name(s): Thrombin receptor
Expression Region	46-432
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.