



Recombinant Human Galanin receptor type 2 (GALR2)

Product Code	CSB-CF009224HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O43603
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MNVSGCPGAGNASQAGGGGGWHPEAVIVPLLALIFLVGTVGNTLVLAVLLRG GQAVSTT NLFILNLGVADLCFILCCVFPFQATIYTLDGWVFGSLLCKAVHFLIFLTMHASSFTL AAVS LDRYLAIKYPLHSRELRTPRNALAAIGLIWGLSLLFSGPYLSYYRQSQLANLTV HPAWS APRRRAMDICTFVFSYLLPVLVGLTYARTLRYLWRAVDPVAAGSGARRAKRK VTRMILI VAALFCLCWMPHHALILCVWFGQFPLTRATYALRILSHLVSYANSCVNPIVYAL VSKHFR KGFRTICAGLLGRAPGRASGRVCAAARGTHSGSVLERESSDLLHMSEAAGAL RPCPGASQ PCILEPCPGPSWQGPAGDSILTVDVA
Source	in vitro E.coli expression system
Target Names	GALR2
Protein Names	Recommended name: Galanin receptor type 2 Short name= GAL2-R Short name= GALR-2
Expression Region	1-387
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.