



Recombinant Macaca mulatta Urotensin-2 receptor (UTS2R)

Product Code	CSB-CF810228MOW
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q8HYC3
Product Type	Transmembrane Protein
Immunogen Species	Macaca mulatta (Rhesus macaque)
Sequence	MALTPESPSSFPGLAAIGSSVPEPPGSPNATLNSSWASPTEPSSLEDLVATGAI GTLLSA MGVVGVVGNAYTLVVTCRSLRAVASMYYVNLALADLLYLLSIPFIVATYITKE WHFGD VGCRVLFSLDFLTMHASIFTLTMSSERYAAVLRPLDTVQRPKG YR KLLALGT WLLALLL TLPVMLAMRLVRRGPKSLCLPAWG PRAH RAYLTLLFATS IAGPGLLIGLLYARL ARAYRR SQRASFKRARRPGARALRLVLGIVLLFWACFLPFWLWQLLAQYREAPLAPRTA RIVNYLT TCLTYGNSCANPFLYTLTRNYRDHLRGRVRS PGSGGVRGPVPSLQPRARFQ RGSGRSLS SCSPQPTESLVLAAAAPAGPALESPGDPA
Source	in vitro E.coli expression system
Target Names	UTS2R
Protein Names	Recommended name: Urotensin-2 receptor Short name= UR-2-R Alternative name(s): Urotensin II receptor Short name= UR-II-R
Expression Region	1-389
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.