



Recombinant *Drosophila melanogaster* Putative odorant receptor 85c (Or85c)

Product Code	CSB-CF893759DLU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q9VHQ6
Product Type	Transmembrane Protein
Immunogen Species	<i>Drosophila melanogaster</i> (Fruit fly)
Sequence	<p>MKFMKYAVFFYTSVGIPTYTIDSRSKKASLWSHLLFWANVINLSVIVFGEILYLG VAYSD GKFIDAVTVLSYIGFVIVGMSKMFFIWWKKTDLSDLVKELEHIYPNGKAEEMY RLDRYL RSCSRISITYALLYSVLIWTFNLF SIMQFLVYEKLLKIRVVGQTLPYLMYFPWNW HENWT YYVLLFCQNFAGHTSASGQISTDLLCAVATQVVMHFDYLARVVEKQVLDRDW SENSRFL AKTVQYHQIRLRLMDVLNDIFGIPLLLNFMVSTFVICFVGFQMTVGVPPDIMIKLF LFLF SSLSQVYLICHYGQLIADASSLSISAYKQNWQNADIRYRRALVFFIARPQRTTY LKATI FMNITRATMTDLLQVSYKFFALLRTMYIK</p>
Source	in vitro <i>E.coli</i> expression system
Target Names	Or85c
Protein Names	Recommended name: Putative odorant receptor 85c
Expression Region	1-389
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>