



Recombinant Macaca mulatta Oxytocin receptor (OXTR)

Product Code	CSB-CF017316MOW
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P56494
Product Type	Transmembrane Protein
Immunogen Species	Macaca mulatta (Rhesus macaque)
Sequence	MEGELAAWSTEAVNSSAAPPGAEGNCTAGPPRRNEALARVEVAVLCLILFLA LSGNACV LLALRTRRHKHSRLFFFMKHLADIADLVVAVFQVLPQLLWDITFRFYGPDLLCRLV KYLQV VGMFASTYLLLLMSLDRCLAICQPLRSLRRRTDRLAVLATWLGCLVASAPQVHI FSLREV ADGVFDCWAVFIQPWGPWAYITWITLAVYIVPVIVLAACYGLISFKIWQNLRLKT AAAAA AEAPEGAAAGDGGRMALARVSSVKLISKAKIRTVKMTFIVLAFIVCWTPFFVQ MWSVW DANAPKEASAFIIVMLLASLNCCNPWIYMLFTGHLFHELVQRFLCCSASYLKG NRLGET STSKKSNSSSFVLSHRSSSQRSQSPSTA
Source	in vitro E.coli expression system
Target Names	OXTR
Protein Names	Recommended name: Oxytocin receptor Short name= OT-R
Expression Region	1-389
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.