



Recombinant Sheep Oxytocin receptor (OXTR)

Product Code	CSB-CF017316SH
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q28756
Product Type	Transmembrane Protein
Immunogen Species	Ovis aries (Sheep)
Sequence	MEGAFANWSAEAVNGSAAPPGTEGNRTAGPPQRNEALARVEVAVLSLILFL ALSGNACV LLALRTTRHKHSRLFFFMKHLSIADLAVAVFQVLPQLLWDITFRFYGPDLLCRLV KYLQV VGMFASTYLLLLMSLDRCLAICQPLRSLRRRTDRLAVLATWLGCLVASAPQVHI FSLREV ADGVFDCWAVFIQPWGPWAYITWITLAVYIVPVIVLAACYGLISFKIWQNLRLKT EAAAA EASAGAEGAAADCAGRAALARVSNVKLISKAKIRTVKMTFIVVLAFIVCWTPFFF KQMWS VWDADAPKEASAFIIAMLLASLNCCNPWIYMLFTGHLFQDLVQRFLCCSFRL KGSQLG ETSVTKKIHSYTFVLSRHSSSQRSCSQPSTV
Source	in vitro E.coli expression system
Target Names	OXTR
Protein Names	Recommended name: Oxytocin receptor Short name= OT-R
Expression Region	1-391
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.