



Recombinant Mouse Somatostatin receptor type 1 (Sstr1)

Product Code	CSB-CF022724MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P30873
Product Type	Transmembrane Protein
Immunogen Species	Mus musculus (Mouse)
Sequence	MFPNGTASSPSSSPSPSPGSCGEGACSRGPGSGAADGMEEPGRNASQNGT LSEGQGSAIL ISFIYSVVCLVGLCGNSMVIYVILRYAKMKTATNIYILNLAIADPELLMLSVPLVTS TLL RHWPFGALLCRLVLSVDAVNMFTSIYCLTVLSVDTRYVAVVHPIKAARYRRPTVA KVVNLG VWVLSLLVILPIVVFSTRTAANS DGT VACNMLMPEPAQRWLVGFLYTFMLMGFLL PVGAIC LCYVLIIAKMRMVALKAGWQQQRKRSEKITLMVMMVVMVFVICWMPFYVQLV NVFAEQD DATVSQSLSVILGYANSCANPILYGFLSDNFKRSFQRILCLSWMDNAAEEPVDYY ATALKS RAYSVEDFQPENLES GG VFRNGTCASRISTL
Source	in vitro E.coli expression system
Target Names	Sstr1
Protein Names	Recommended name: Somatostatin receptor type 1 Short name= SS-1-R Short name= SS1-R Short name= SS1R Alternative name(s): SRIF-2
Expression Region	1-391
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.