



Recombinant Chicken Corticotropin-releasing factor receptor 1 (CRHR1)

Product Code	CSB-CF841911CH
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q90812
Product Type	Transmembrane Protein
Immunogen Species	Gallus gallus (Chicken)
Sequence	SIQEYCESLLPTTNHTGPQCNASVDLIGTCWPRSAVGQLVARPCPEYFYGV RYNTTNG YRECLANGSWAARVNYSQCQEILSEEKRSKLHYHIAVIINYLGHCVSLGTLLVA FVLFMR LRSIRCLRNIHWNLITAFILRNATWFVQVQLTMNPEVHESNVVWCRLVTAAANY FHVTNF FWMFGEGCYLHTAIVLTYSTDKLRKWMFICIGWCIPFPPIIWAWAIGKLYDNEKC WFGKR AGVYTDYIYQGPMILVLLINFILFNIVRILMTKLRASSTSETIQYRKAVKATLVLLS LL GITYMLFFVNPGEDEISRIVFIYFNSFLESFQGFFVSVFYCFLNSEVRS AVRKRW HRWQD KHSIRARVARAMSIPTSPTRVVSFHSIKQSSAV
Source	in vitro E.coli expression system
Target Names	CRHR1
Protein Names	Recommended name: Corticotropin-releasing factor receptor 1 Short name= CRF-R-1 Short name= CRF-R1 Short name= CRFR-1 Alternative name(s): Corticotropin-releasing hormone receptor 1 Short name= CRH-R-1 Short name=
Expression Region	29-420
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.