



Recombinant Human Prokineticin receptor 1 (PROKR1)

Product Code	CSB-CF018748HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q8TCW9
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	METTMGMFMDNATNTSTSFSLSVLNPHGAHATSFPFNFSYSDYDMPLEDEDV TNSRTFFA AKIVIGMALVGIMLVCGIGNFIFIAALVRYKCLRNLNLLIANLAISDFLVAIVCCPF EM DYYVVRQLSWEHGHVLCSTSVNYLRTVSLYVSTNALLAIAIDRYLAIVHPLRPRM KCQTAT GLIALVWTVSILIAIPSAYFTTETVLVIVKSQEKIFCGQIWPVDQQLYYKSYFLFIF GIE FVGPVVTMTLCYARISRELWFKAVPGFQTEQIRKRLRCRRKTVLVLMCILTAYV LCWAPF YGFTIVRDFFPVVFVKEKHLYLTA FYIVECIAMSNMINTLCFVTVKNDTVKYFKKI MLLH WKASYNGGKSSADLDLKTIGMPATEEVDCIRLK
Source	in vitro E.coli expression system
Target Names	PROKR1
Protein Names	Recommended name: Prokineticin receptor 1 Short name= PK-R1 Alternative name(s): G-protein coupled receptor 73 G-protein coupled receptor ZAQ GPR73a
Expression Region	1-393
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.