



Recombinant Mouse Prokineticin receptor 1 (Prokr1)

Product Code	CSB-CF018748MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q9JKL1
Product Type	Transmembrane Protein
Immunogen Species	Mus musculus (Mouse)
Sequence	METTVGALGENTTDTFTDFFSALDGHEAQTGSLPFTFSYGDYDMPLDEEEDV TNSRTFFA AKIVIGMALVGIMLVCGIGNFIFITALARYKKLRNLTNLLIANLAISDFLVAIVCCPF EM DYYVVRQLSWEHGHVLCASVNYLRTVSLYVSTNALLAIAIDRYLAIVHPLRPRM KCQTAA GLIFLVWSVSILIAIPAAYFTTETVLVIVERQEKIFCGQIWPVDQQFYYSYFLLVF GLE FVGPPVAMTLCYARVSRELWFKAVPGFQTEQIRRLRCRRRTVLGLVCVLSA YVLCWAPF YGFTIVRDFFPSVFVKEKHLYTAFYVVECIAMSNMINTLCFVTVRNNTSKYLK RILRLQ WRASPSGSKASADLDLRTTGIPATEEVDCIRLK
Source	in vitro E.coli expression system
Target Names	Prokr1
Protein Names	Recommended name: Prokineticin receptor 1 Short name= PK-R1 Alternative name(s): G-protein coupled receptor 73
Expression Region	1-393
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.