



Recombinant *Talpa europaea* Alpha-2B adrenergic receptor (ADRA2B)

Product Code	CSB-CF001389TLA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O19091
Product Type	Transmembrane Protein
Immunogen Species	<i>Talpa europaea</i> (European mole)
Sequence	AIAAIIIFLILFTIFGNALVILAVLTSRSLRAPQNLFLVSLAAADILVATLIIPFSLANE LLGYWYFRRMWCKVYLALDVLFCSTSSIVHLCAISLDRYWAVSRALEYNSKRTP RRIKCII LMVWLIAAVISLPSLVYKGDQGPQPSGAPQCNLNQETWYILASSIGSFFAPCLI MILVYL RIYLIAKRSHCRGPRAKGAPGKSKFKQSRQVPGGTLASAKVPNLA AHLVAAGE TNGRSKP TGEKEMGETPEDSGTSTLPPSWPALPNSGQDQKEGICEASLEEEAE EEEEEGE E EEEEECE PQALPASPASACSPPLQQPQGSQVLATLRGQVLLGRAAGAASGQWRRRTQ LSREKRFTF VLAVVIGVFVLCWFPFFFSYSLGAICPQQCKVPHDLF
Source	in vitro E.coli expression system
Target Names	ADRA2B
Protein Names	Recommended name: Alpha-2B adrenergic receptor Alternative name(s): Alpha-2B adrenoreceptor Short name= Alpha-2B adrenoceptor Short name= Alpha-2BAR
Expression Region	1-397
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.