



Recombinant Drosophila melanogaster Odorant receptor 22b (Or22b)

Product Code	CSB-CF305927DLU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P81910
Product Type	Transmembrane Protein
Immunogen Species	Drosophila melanogaster (Fruit fly)
Sequence	MLSQFFPHIKEKPLSERVKSRDAFVYLDRVMWSFGWTVPENKRWDLHYKLW STFVTLIF ILLPISVSVEYIQRFKTFSAGEFLSSIQIGVNMYGSSFKSYLTMMGYKKRQEAKM SLDEL DKRCVCDEERTIVHRHVALGNFCYIFYHIA YTSFLISNFLSFIMKRIHAWRMYFP YVDPE KQFYISSIAEVILRGWAVFMDLCTDVCPLISMVIARCHITLLKQRLRNL RSEPGR TEDEY LKELADCVRDHR LILDYVDALRSVFSGTIFVQFLLIGIVLGLSMINIMFFSTLSTGV AVV LFMSCVSMQTFPFCYLCNMIMDDCQEMADSLFQSDWTSADRRYKSTLVYFLH NLQQPIIL TAGGVFPISMQTNLNMVKLAFTVVTIVKQFNLAEKFQ
Source	in vitro E.coli expression system
Target Names	Or22b
Protein Names	Recommended name: Odorant receptor 22b
Expression Region	1-397
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.