



Recombinant *Drosophila melanogaster* Putative gustatory receptor 59e (Gr59e)

Product Code	CSB-CF895227DLU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q9W1N6
Product Type	Transmembrane Protein
Immunogen Species	<i>Drosophila melanogaster</i> (Fruit fly)
Sequence	<p>MDSSYWENLLLITINRFLGVYPSGRVGVLRWLHTLWSLFLLMYIWTGSIVKCLE FTVEIPT IEKLLYLMEFPGNMATIAILVYYAVLNRPLAHGAELQIERIITGLKGGKAKRLVYKR HGQR TLHLMATTLVFHGLCVLVDVVNYDFEFWTTWSSNSVYNLPGLMMSLGVLYA QPVHFLWL VMDQMRMCLKELKLLQRPPQGSTKLDACYESAFVLDAGGGSSALMIEEMRY TCNLIEQV HSQFLFRFGLYLVNLLNSLVSICVELYLIFNFFETPLWEESVLLVYRLLWLAMH GGRIW FILSVNEQILEQKCNLCQLLNELEVCSRLQRTINRFLQLQRSIDQPLEACGIVT LDTR SLGGFIGVLMIAIVFLIQIGLGNKSLMGVALNRSNWVYV</p>
Source	in vitro E.coli expression system
Target Names	Gr59e
Protein Names	Recommended name: Putative gustatory receptor 59e
Expression Region	1-399
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>