



Recombinant Sheep Bombesin receptor subtype-3 (BRS3)

Product Code	CSB-CF002820SH
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O97967
Product Type	Transmembrane Protein
Immunogen Species	Ovis aries (Sheep)
Sequence	MSQRQPQSPNQTLISTTNDTESSSSVVPNDSTNKRRTGDNSPGIEALCAIYITY AVIISV GILGNAILIKVFFKTKSMQTPNIFITSLAFGDLLLLLTCVPVDVTHYLAEGWLF RIGC KVLFSIRLTSVGVSVFTLTILSADRYKAVVKPLERQPPNAILKTCAKAGCIWIMS MIAL PEAIFSNVYTFQDPDKNVTFKACASYPVSERLLQEIHSLCFLVFYIIPLSISVYY SLI ARTLYKSTLNIPTEEQRHARKQIESRKRIAKTVLVLVALFALCWLPNHLLYLYRS FTSQT YMDSSTVHLFVTIISRILAFSNSCVNPFALYWLSNTFQQHFKAQLFCCKAGRPD PTAANT PLDNLAVMGRVPGAASTQMSEISVSPFTGCSVKKEDDRV
Source	in vitro E.coli expression system
Target Names	BRS3
Protein Names	Recommended name: Bombesin receptor subtype-3 Short name= BRS-3
Expression Region	1-399
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.