



Recombinant Pig Mu-type opioid receptor (OPRM1)

Product Code	CSB-CF842513PI
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q95247
Product Type	Transmembrane Protein
Immunogen Species	Sus scrofa (Pig)
Sequence	MDSSADPRNASNCTDPFSPSSMCSPVPSPSSWVNFHLEGNLSDPCIRNRTE LGGSDSLC PPTGSPSMVTAITIMALYSIVCVVGLFGNFLVMYVIVRYTKMKTATNIYIFNLALA DALA TSTLPFQSVNYLMGTWPFGTILCKIVISIDYYNMFTSIFTLCTMSVDRIYAVCHPW KALD FRTPRNAKIINVCNWILSSAIGLPVMFMATTKYRNGSIDCALTFSHPTWYWENL LKICVF IFAFIMPVLIITVCYGLMILRLKSVRMLSGSKEKDRNLRRITRMVLVVAVFIVCW TPIH IYVIIKALITIPETTFQTVSWHFCIALGYTNSCLNPVLYAFLDENFKRCRFREFCIPT SST IEQQNSARIRQNTRDHPSTANTVDRTNHQLENLEAETAPLP
Source	in vitro E.coli expression system
Target Names	OPRM1
Protein Names	Recommended name: Mu-type opioid receptor Short name= M-OR-1 Short name= MOR-1
Expression Region	1-401
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.