



Recombinant Bovine Mu-type opioid receptor (OPRM1)

Product Code	CSB-CF016361BO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P79350
Product Type	Transmembrane Protein
Immunogen Species	Bos taurus (Bovine)
Sequence	MDSGAVPTNASNCTDPFTHPSSCSPAPSPSSWVNFHLEGNLSDPCGPNRT ELGGSDRLC PSAGSPSMITAIIMALYSIVCVVGLFGNFLVMYVIVRYTKMKTATNIYIFNLALAD ALA TSTLPFQSVNYLMGTWPFGTILCKIVISIDYYNMFTSIFTLCTMSVDRIYAVCHPV KALD LRTPRNAKIINICNWILSSAIGLPVMFMATTKYRQGSIDCTLTFSHPTWYWENLL KICVF IFAFIMPILITVCYGLMILRLKSVRMLSGSKEKDRNLRRITRMVLVVAVFIVCWT PIH IYVIIKALITIPETTFQTVSWHFCIALGYTNSCLNPVLYAFLDENFKRCRFREFCIPT SST IEQQNSTRIRQNTRDHPSTANTVDRTNHQLENLEAETPLP
Source	in vitro E.coli expression system
Target Names	OPRM1
Protein Names	Recommended name: Mu-type opioid receptor Short name= M-OR-1 Short name= MOR-1
Expression Region	1-401
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.