



Recombinant Pig Growth hormone-releasing hormone receptor (GHRHR)

| | |
|--------------------------|--|
| Product Code | CSB-CF009413PI |
| Storage | Store at -20°C, for extended storage, conserve at -20°C or -80°C. |
| Uniprot No. | P34999 |
| Product Type | Transmembrane Protein |
| Immunogen Species | Sus scrofa (Pig) |
| Sequence | HVHPECDFITQLREDERTCLQAADRMANSSSGCPRTWDGLLCWPTAGPGEW VTLPCPAFF SHFSSEPGALKRDCTTTGWSEPFPPYPEACPVPLELLTDEKSYFSTVRIVYTT GHSVSAV ALFVAIAILVALRRLHCPRNYIHSQLFATFILKAGAVFLKDAALFHSENTDHCSFS TVLC KVS VATSHFATMTNFSWLLAEAVYLTCLLASTSPSTRRAFWWLVLAGWGLPLL FTGTWVG CKLAFEDVACWDLDDSSPYWWIKGPIVLSVGVNFGLFLNIIRILLRKLEPAQGS LHTQP QYWRLSKSTLLLIPLFGIHYVIFNFLPDSAGLGIRLPLELGLGSFQGFIVAILYCFL NQE VRTEISRRWHGHDPELLPAWRTHAKWAKPSRSRAKVLTTVC |
| Source | in vitro E.coli expression system |
| Target Names | GHRHR |
| Protein Names | Recommended name: Growth hormone-releasing hormone receptor Short name= GHRH receptor Alternative name(s): Growth hormone-releasing factor receptor Short name= GRF receptor Short name= GRFR |
| Expression Region | 23-423 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | Full Length of Mature Protein |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |