



Recombinant Guinea pig Substance-K receptor (TACR2)

Product Code	CSB-CF737329GU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q64077
Product Type	Transmembrane Protein
Immunogen Species	Cavia porcellus (Guinea pig)
Sequence	MGACVIVTNTNISSGLESNTTGITAFSMPTWQLALWATAYLALVLVAVTGNATV TWILA HQRMRVTNYFIVNLALADLCMAAFNAAFNFVYASHNIWYFGRAFCYFQNLFP TAMFVS IYSMTAIAIDRYMAIVHPFQPRLSAPSTKAVIGGIWLVALALAFPQCFYSTITEDE GATK CVVAWPEDSRDKSLLLYHLVVIVLIYLLPLTVMFVAYSIIGLTLWRRRAVPRHQA GANLR HLQAKKKFVKTMVLVVVTFaicwLpYHLYFILGSFQEDIYCHKFIQQVYLALFWL AMSST MYNPIIYCCLNRRFRSGFRLAFRCCPWVTPTEEDKLELTHTPSFSLRVNRCHT KEILFMA GDTV PSEATNGQAGGPQDRESVELSSLPGCRAGPSILAKASS
Source	in vitro E.coli expression system
Target Names	TACR2
Protein Names	Recommended name: Substance-K receptor Short name= SKR Alternative name(s): NK-2 receptor Short name= NK-2R Neurokinin A receptor Tachykinin receptor 2
Expression Region	1-402
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.