



Recombinant Rat Prostaglandin D2 receptor 2 (Ptgdr2)

Product Code	CSB-CF744318RA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q6XKD3
Product Type	Transmembrane Protein
Immunogen Species	Rattus norvegicus (Rat)
Sequence	<p>MANITLKPLCPLLEEMVQLPNHSNSSLRYIDHVSVLLHGLASLLGLVENGLILFV VGCRM RQTVVTTWVLHLALSDLLAAASLPFFTYFLAVGHSWELGTTFCCLHSSVFFLN MFASGFL LSAISLDRCLQVVRPVWAQNHRTVAAAHRVCLMLWALAVLNTVPYFVFRDITP RRDGRIM CYYNMLLLNPGSDRDTTCDYRQKALAVSKFLLAFMVPLAIIASSHVAVSLQLHH RGRQRT GRFVRLVAAIVVAFILCWGPYHIFSLLEARAHSVTTLRQLASRGLPFVTSLAFFN SVVNP LLYVLTCPDMLHKLRRSLLTVLESVLVEDSDLSTGPGKRCRRRHRRRASSTTT PASTLLL ADRFPQLRPARLIGWMRRGSAELPRRVREQSQEKGQGLSCTLD</p>
Source	in vitro E.coli expression system
Target Names	Ptgdr2
Protein Names	Recommended name: Prostaglandin D2 receptor 2 Alternative name(s): G protein-coupled receptor 44
Expression Region	1-403
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.