



Recombinant Drosophila melanogaster Odorant receptor 43b (Or43b)

Product Code	CSB-CF305793DLU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P81918
Product Type	Transmembrane Protein
Immunogen Species	Drosophila melanogaster (Fruit fly)
Sequence	MFGHFKLVYPAPISEPIQSRDSNAYMMETLRNSGLNLKNDFGIGRKIWRVFSFT YNMVIL PVSFPINIVIHIAEFPPPELLLQSLQLCLNTWCFALKFFTLIVYTHRLELANKHFDE LDKY CVKPAEKRKVRDMVATITRLYLTFVVVYVLYATSTLLDGLLHHRVPYNTYYPFIN WRVDR TQMYIQSFLEYFTVGYAIYVATATDSYPVIYVAALRTHILLKDRIIYLGDPSSNEG SSDP SYMFKSLVDCIKAHRTMLNFCDAIQPIISGTIFAQFIICGSILGIIMINMVLFAQDQST RF GIVYVMAVLLQTFPLCFYCNAIVDDCKELAHALFHSAWWWQDKRYQRTVIQFL QKLQQP MTFTAMNIFNINLATNINVAKFAFTVYAIASGMNLDQKLSIKE
Source	in vitro E.coli expression system
Target Names	Or43b
Protein Names	Recommended name: Odorant receptor 43b
Expression Region	1-403
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.