



Recombinant Protein sbmA (sbmA)

Product Code	CSB-CF360344EOD
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P0AFY7
Product Type	Transmembrane Protein
Immunogen Species	Escherichia coli O157:H7
Sequence	MFKSFFPKPGTFFLSAFVWALIAVIFWQAGGGDWVARITGASGQIPISAARFW SLDFLIF YAYYIVCVGLFALFWFIYSPHRWQYWSILGTALIIFVTWFLVEVGVAVNAWYAP FYDLIQ TALSSPHKVTIEQFYREVGVLGIALIAVVISVLNNFFVSHYVFRWRTAMNEY Y MANWQQ LRHIEGAAQRVQEDTMRFASTLENMGVSFINAIMTLIAFLPVLVTL SAHVPELP II GHIP YGLVIAAIVWSLMGTGLLAVVGIKLPGLEFKNQRVEAAYRKELVYGEDD ATRAT PPTVRE LFSAVRKNYFRLYFHMYFNIARILYLQVDNVFGLFLLFPSIVAGTITLGLMT QIT NVFG QVRGAFQYLINSWTTLVELMSIYKRLRSFEHELDGDKIQEVTH TLS
Source	in vitro E.coli expression system
Target Names	sbmA
Protein Names	Recommended name: Protein sbmA
Expression Region	1-406
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.