



Recombinant Dog Probable G-protein coupled receptor 83 (GPR83)

Product Code	CSB-CF879996DO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q9TTQ9
Product Type	Transmembrane Protein
Immunogen Species	Canis lupus familiaris (Dog) (Canis familiaris)
Sequence	AERPEGRADEPGLAALAGPNASHFFWSNYSFSDWQNFVGRRRRYGAESQNP TVKALLVVA YSFIIVFSLFGNVLVCHVIFKNQRMRSATSLFIVNLAVADILITLLNTPFTLVRFVN STW VFGKGMCHVSRFAQYCSLHVSALTLTAIIVDRHQVIMHPLKPRISITKGVIIYITVI WTMA TFFSLPHAICQKLFTFKYSEDIVRSLCLPDFPEPADLFWKYLDLATFILLYILPLLI SV AYARVAKKLWLCNTIGDVTTEQYLALRRKKKTIKMLMLVVVLFALCWFPPLNCY VLLLSS KVIHTNNALYFAFHWFAMSSTCYNPFIYCWLNENFRIELKALLSMCQRLPKPQE ERPPSP VPSFRVAWTEKSNRRVPPANLLLSSHLSHGKTDLSSVEPIVAMS
Source	in vitro E.coli expression system
Target Names	GPR83
Protein Names	Recommended name: Probable G-protein coupled receptor 83 Alternative name(s): Glucocorticoid-induced receptor
Expression Region	17-422
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.