



Recombinant Drosophila melanogaster Putative odorant receptor 10a (Or10a)

Product Code	CSB-CF892919DLU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q9VYZ1
Product Type	Transmembrane Protein
Immunogen Species	Drosophila melanogaster (Fruit fly)
Sequence	MSEWLRFLKRDQQLDVYFFAVPRLSLDIMGYWPGKTGDTWPWRSLIHFAILAI GVATELH AGMCFDRQQITLALETLCPAGTSAVTLKMFMLRFRQDLSIMWNRLRGLLF DPNWERP EQRDIRLKHSAMAARINFWPLSAGFFTCTTYNLKPILIAMILYLQNRVEDFWFT PFNMT MPKVLLNYPFFPLTYIFIAYTYGYVTIFMFGGCDGFYFEFCAHLSALFEVLQAEIE SMFRP YTDHLELSPVQLYILEQKMRSVIIRHNAIDLTRFFRDRYTIITLAHFVSAAMVIGF SMV NLLTLGNNGLGAMLYVAYTVAALSQLLVYCYGGTLVAESSTGLCRAMFSCPW QLFKPKQR RLVQLLILRSQRPVSMVAVPFFSPSLATFAAILQTSGSIIALVKSFQ
Source	in vitro E.coli expression system
Target Names	Or10a
Protein Names	Recommended name: Putative odorant receptor 10a
Expression Region	1-406
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.