



Recombinant Drosophila melanogaster Putative odorant receptor 67a (Or67a)

Product Code	CSB-CF895188DLU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q9VT08
Product Type	Transmembrane Protein
Immunogen Species	Drosophila melanogaster (Fruit fly)
Sequence	MDNVAEMPEEKYVEVDDFLRLAVKFYNTLGIDPYETGRKRTIWFQIYFALNMF NMVFSFY AEVATLVDRLRDNENFLESCILLSYVSFVVMGLSKIGAVMKKKPKMTALVRQLE TCFPSP SAKVQEEYAVKSWLKRCHIYTKGFGGLFMIMYFAHALIPLFIYFIQRVLLHYPDA KQIMP FYQLEPWEFRDSWLFYPSYFHQSSAGYTATCGSIAGDLMIFAVVLQVIMHYER LAKVLRE FKIQAHNAPNGAKEDIRKLQSLVANHIDILRLTDLMNEVFGIPLLLNFIASALLVCL VGV QLTIALSPEYFCKQMLFLISVLLLEVYLLCSFSQRLIDASENVGHAAYDMDWLGS DKRFKK ILIFISMRSQKPVCLKATVVLDLSMPTMSIFLGMYSYKFFCAVRTMYQ
Source	in vitro E.coli expression system
Target Names	Or67a
Protein Names	Recommended name: Putative odorant receptor 67a
Expression Region	1-407
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.