



Recombinant Rat Substance-P receptor (Tacr1)

Product Code	CSB-CF023068RA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P14600
Product Type	Transmembrane Protein
Immunogen Species	Rattus norvegicus (Rat)
Sequence	MDNVLPMDSDLFPNISTNTSESNQFVQPTWQIVLWAAAYTVIVVTSVVGNNVVI WIILAH KRMRTVTNYFLVNLAFAEACMAAFNTVVNFTYAVHNVWYYGLFYCKFHNFFPI AALFASI YSMTAVAFDRYMAIIHPLQPRLSATATKVVFVIWVLALLLAFPQGYSTTETMP SRVVC MIEWPEHPNRTYEKAYHICVTVLIYFLPLLIGYAYTVVGITLWASEIPGDSSDR YHEQV SAKRKVVKMMIVVCTFAICWLPFHVFFLLPYINPDLYLKKFIQQVYLASMWLA MSSTMY NPIIYCCLNDRFRLGFKHAFRCCPFISAGDYEGLEMKSTRYLQTQSSVYKVSRL ETTIST VVGAHEEEPEEGPKATPSSLDLTSNGSSRSNSKTMTESSSFYSNMLA
Source	in vitro E.coli expression system
Target Names	Tacr1
Protein Names	Recommended name: Substance-P receptor Short name= SPR Alternative name(s): NK-1 receptor Short name= NK-1R Tachykinin receptor 1
Expression Region	1-407
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.