



Recombinant Danio rerio Alpha-2Da adrenergic receptor (adra2da)

Product Code	CSB-CF811652DIL
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q8JG70
Product Type	Transmembrane Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Sequence	MSVTPTANSTEEAANITASPRLWPYTEPASAIILLVSLIILLTIVGNVLVIVAVLTS RA LRAPQNLFLVSLACADILVATLVIPFLANEIMGYWYFGSTWCAFYALDLVLFCT SSIVH LCAISLDYWSVTKAVRYNLKRTPRRIKCMIAVVWLISAVISFPPLIMTKHDEKE CLIND ETWYILSSCAVSFFAPGLIMITVYCKIYRVAKQRSSTVFVAKNGLERQPSQSET CFVRKD KFEKESPSSNSSESAQRQEELDDIDLEESAASDSRARGSRFSKRRRVEGERR GPQRSCRV SWAAHQEPGSRQQQLASKSKVAQMREKRFTFVLAVVMGVFVLCWFPPFFFTY SLQAVCGER CGPPEALFKLFFWIGYCNSSVNPIIYTIFNRDFRKAFKKVVCWSAQRT
Source	in vitro E.coli expression system
Target Names	adra2da
Protein Names	Recommended name: Alpha-2Da adrenergic receptor Alternative name(s): Alpha-2Da adrenoceptor Short name= Alpha(2Da)AR Alpha-2Da adrenoreceptor
Expression Region	1-408
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.