



Recombinant *Drosophila melanogaster* Putative odorant receptor 65c (Or65c)

Product Code	CSB-CF306593DLU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P82984
Product Type	Transmembrane Protein
Immunogen Species	<i>Drosophila melanogaster</i> (Fruit fly)
Sequence	<pre>MDIRGNVHRFVKFYIDGWKHFRDPTMESSYSAVYYWREQMKAMFLYTTSKER QMPYRSSW HTLVIIQATVCFLTMCYGVTESLGDKVQMGRDIAFIIGFFYIAFKIYYFQWYGDEL DEVV EALETFHPWAQKGPQAVDYRTAKRWYFTLAFFLASSWLFLCIFILLITSPW VHQQIL PLHAAFPPQWHEKSIHPISHAFIYLFQTNVVMYFLTWLVCIEGLSVSIYVEITFAI EVLC LELRHLHQRCHGYEQLRLETNRLVQFHQKIVHILDHTNKVFHGTLMQMGVNF FLVSLSV LEAMEARKDPKVVAQFAVLMMLLALGHLSMWSYFGDLLSQKSLTISEAAYEAYD PIKGSKD VYRDLCLIIRRGQEPLIMRASPFPFNFNYSAILNQCYGILTFLKTLTD</pre>
Source	in vitro <i>E.coli</i> expression system
Target Names	Or65c
Protein Names	Recommended name: Putative odorant receptor 65c
Expression Region	1-410
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.