



Recombinant *Drosophila melanogaster* Putative odorant receptor 59c (Or59c)

Product Code	CSB-CF896182DLU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q9W1P7
Product Type	Transmembrane Protein
Immunogen Species	<i>Drosophila melanogaster</i> (Fruit fly)
Sequence	MTKFFFKRLQTAPLDQEVSSLDASDYYYRIAFFLGWTPPKGALLRWIYSLWTLT TMWLG VYLPLGLSLTYVKHFDRFTPTEFLTSLQVDINCIGNVIKSCVTYSQMWRFRMN ELISSL DKRCVTTTQRRIFHKMVARVNLIVILFLSTYLGFCLTLFTSVFAGKAPWQLYNP LVDWR KQHWQLWIASILEYCVVSIGTMQELMSDTYAIVFISLFRCHLAILRDRIANLRQD PKLSE MEHYEQMVACIQDHRITIIQCSQIIRPILSITIFAQFMLVGIDLGLAAISILFFPNTIW TI MANVSFIVAICTESFPCCMLCEHLIEDSVHVSNALFHSNWITADRSYKSAVLYFL HRAQQ PIQFTAGSIFPISVQSNIAVAKFAFTIITIVNQMNLGKFFSDRSNGDINP
Source	in vitro <i>E.coli</i> expression system
Target Names	Or59c
Protein Names	Recommended name: Putative odorant receptor 59c
Expression Region	1-411
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.