



Recombinant Dog Adenosine receptor A2a (ADORA2A)

Product Code	CSB-CF001376DO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P11617
Product Type	Transmembrane Protein
Immunogen Species	Canis lupus familiaris (Dog) (Canis familiaris)
Sequence	MSTMGSWVYITVELAIAVLAILGNVLCWAVWLNSNLQNVNTNYFVVSLLAAADIA VGVLAI PFAITISTGFCAACHNCLFFACFVLVLTQSSIFSLLAIAIDRYIAIRIPLRYNGLVTG TR AKGIIAVCWVLSFAIGLTPMLGWNNCSQPKEGRNYSQGC GEGQVACLFEDVV PMNYMVYY NFFAFVLVPLLLMLGVYLRIFLAARRQLKQMESQPLPGERARSTLQKEVHAAK SLAIIVG LFALCWLPLHIINCFTFFCPEC SHAPLWLMYLTIVLSHTNSVVPFIYAYRIREFR QTFR KIIRSHVLR RREPFKAGGTSARALAAHGSDGEQISLRLNGHPPGVWANGSAPH PERRPNG YTLGLVSGGIAPESHGDMGLPDVELLSHELK GACPE SPGLEGPLAQDGAGVS
Source	in vitro E.coli expression system
Target Names	ADORA2A
Protein Names	Recommended name: Adenosine receptor A2a
Expression Region	1-412
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.