



Recombinant Mouse Cytochrome c oxidase assembly protein COX15 homolog (Cox15)

Product Code	CSB-CF804356MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q8BJ03
Product Type	Transmembrane Protein
Immunogen Species	Mus musculus (Mouse)
Sequence	<p>MQRLLPSLRALTGSRNVGLLVPRVASRTQCGSSCGSQRPLRPGQYSTITEVA LQSGKGT VSLPSRAAERAVGRWLLVCSGTVAGAVILGGVTRLTESGLSMVDWHLIKEMKP PTSQEEW EAEFQKYQQFPEFKILNHDMTLAEFKFIWYMEYSHRMWGRAVGLAYILPAAFY WRKGWLN RGMKGRVLALCGLVCFQGLLGWYMVKSGLEEKPEFYDIPRVSQYRLAAHLGS ALVLYCAS LWTSLSLLLPQHKLPELRQLLWLRRFAGGTAGLVFLTALSGAFVAGLDAGLVY NSFPKMG ESWIPEDLLTFSPVLKNVFENPTMVQFDHRLLGVTSTAITVLYFLSRRIPRR TKMAA VTLLALAYAQVALGISTLLMYVPTPLAATHQSGSLALLSGALWLMNELRRVPK</p>
Source	in vitro E.coli expression system
Target Names	Cox15
Protein Names	Recommended name: Cytochrome c oxidase assembly protein COX15 homolog
Expression Region	1-413
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.