



Recombinant Rat Probable G-protein coupled receptor 19 (Gpr19)

Product Code	CSB-CF009793RA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P70585
Product Type	Transmembrane Protein
Immunogen Species	Rattus norvegicus (Rat)
Sequence	MGFDHRMETDQPPVVTATLLVPLQNGSCVEAAEALLPHGLMELHEEHSWMS NRTDLQYEL NPGEVATASIFFGALWLFSIFGNSLVCLVIHRSRRTQSTTNYFVVSMACADLLIS VASTP FVVLQFTTGRWTLGSAMCKVVRYPQYLTPGVQIYVLLSICIDRFYTIVYPLSFKV SREKA KRMIAASWILDAAFVTPVFFFYGSNWDSHCNYFLPPSWEGTAYTVIHFLVGFVI PSVLII LFYQKVIKIWRIGTDGRTLRRMTNIVPRTKVKTVMFLLLNLVFLFSWLPFHVA QLWHP HEQDYRKSSLVFTAVTWVSFSSSASKPTLYSIYNANFRRGMKETFCMSSMKC YRSNAYTI TTSSRMAKRNYVGISEIPPVSRITITKDSIYDSFDREAREKKLAWPINSNPPNTFV
Source	in vitro E.coli expression system
Target Names	Gpr19
Protein Names	Recommended name: Probable G-protein coupled receptor 19
Expression Region	1-415
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.