



Recombinant Bovine Endothelin B receptor (EDNRB)

Product Code	CSB-CF007404BO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P28088
Product Type	Transmembrane Protein
Immunogen Species	Bos taurus (Bovine)
Sequence	EEREFPPAGATQPLPGTGEMMETPTETSWPGRSNASDPRSSATPQIPRGGR MAGIPRTP PPCDGPIEKETFKYINTVVSLVFLGIIGNSTLLRIYKKNKCMRNGPNILIASLAL GD LLHIIIDIPINTYKLLAKDWPFVEMCKLVFPIQKASVGITVLSLICALSIDRYRAVA SWS RIKGIGVPKWTAVEIVLIWVSVVLAVPEAVGFDIITSDHIGNKLRICLLHPTQKT AFMQ FYKTAKDWWLFSFYFCLPLAITALFYTLMTCEMLRKKSGMQIALNDHLKQRRE VAKTVFC LVLVFALCWLPLHLSRILKLTLYDQHDP RRCEFLSFLLVLDYIGINMASLNSCINP IALY LVSKRFKNCFKSCLCCWCQSFEKQSLEEKQSCLKFKANDHGYDNFRSSNKY SSS
Source	in vitro E.coli expression system
Target Names	EDNRB
Protein Names	Recommended name: Endothelin B receptor Short name= ET-B Short name= ET-BR Alternative name(s): Endothelin receptor non-selective type
Expression Region	27-441
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.