



Recombinant Pig Endothelin B receptor (EDNRB)

Product Code	CSB-CF007404PI
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P35463
Product Type	Transmembrane Protein
Immunogen Species	Sus scrofa (Pig)
Sequence	EERGFPPAGATPPALRTGEIVAPPTKTFWPRGSNASLPRSSSPQMPKGGRM AGPPARTL TPPPCEGPIEIKDTFKYINTVVSCLVFVLGIIGNSTLLRIIYKNKCMRNGPNILIASL AL GDLLHIIIDIPINVYKLLAEDWPFVEMCKLVFFIQKASVGITVLSLCALSIDRYRA VAS WSRIKGIGVPKWTAVEIVLIWVVSVVLAVPEALGFDMITTDYKGNRLRICLLHPT QKTAF MQFYKTAKDWWLFSFYFCLPLAITAFFYTLMTCEMLRKKSGMQIALNDHLKQR REVAKTV FCLLVLFALCWLPPLHLSRILKLTLYDQNDNSNRCELLSFLLVLDYIGINMASLNSCI NPIA LYLVSKRFKNCFKSCLCCWCQSFEKQSLKFKANDHGYDNFRSSN KYSSS
Source	in vitro E.coli expression system
Target Names	EDNRB
Protein Names	Recommended name: Endothelin B receptor Short name= ET-B Short name= ET-BR Alternative name(s): Endothelin receptor non-selective type
Expression Region	27-443
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.