



Recombinant Human Somatostatin receptor type 3 (SSTR3)

Product Code	CSB-CF022726HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P32745
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	<p>MDMLHPSSVSTTSEPENASSAWPPDATLGNVSAGPSPAGLAVSGVLIPLVYLV VCVVGLL GNSLVIYVVL RHTASPSVTNVYILNLALADELFMLGLPFLAAQNALSYWPFGL MCRLVM AVDGINQFTSIFCLTVMSVD RYLAVVHPTRSARWRTAPVARTVSAAVWVASAV VVLPPVVV FSGVPRGMSTCHMQWPEPAAAWRAGFIIYTAALGFFGPLLVICLCYLLIVVKVR SAGRRV WAPSCQRRRRSERRVTRMVVAVALFVLCWMPFYVLNIVNVVCP LPEEPAFF GLYFLVVA LPYANSCANPILYGFLSYRFKQGFRRVLLRPSRRVRSQEPTVGPPEKTEEEDE EEEDGE SREGGKGKEMNGRVSQITQPGTSGQERPPSRVASKEQQLLPQEASTGEKSS TMRISYL</p>
Source	in vitro E.coli expression system
Target Names	SSTR3
Protein Names	Recommended name: Somatostatin receptor type 3 Short name= SS-3-R Short name= SS3-R Short name= SS3R Alternative name(s): SSR-28
Expression Region	1-418
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.