



Recombinant *Ashbya gossypii* ADIPOR-like receptor IZH3 (IZH3)

Product Code	CSB-CF748429DOT
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q75F81
Product Type	Transmembrane Protein
Immunogen Species	<i>Ashbya gossypii</i> (strain ATCC 10895 / CBS 109.51 / FGSC 9923 / NRRL Y-1056) (Yeast) (<i>Eremothecium gossypii</i>)
Sequence	<p>MSHPNTHMPRTHAVHGRAAPQRRGCRSVEKTDQSGHSRSSLAESAMEQA QLRSRGEAGG GRSVLCASGEPGSAHKWAGAVTCSGAEEQPLHWAESRARGLARHLHYWELP YAWRENYI IYGHRFYHSHRKSLLSVLNAYGWHNETINIWSHLVGA AVLAYLLCWGWPRSDV YRAAQVP RLAKWAIGAFLACGVKCMASVAVWHTFNGTCHLKLRSRVCVDYTGITLLVTA SVVTTVA VTLYGLSRPLMYAYMVASIGLGTAAGVMNWSPHFDRPEARPLRIAVYVGLAAL GLVSFVH VWMQVRWASAHLMAPLVYKSLVWYGIGVVFYATLPERWRSVDVTLDCSSGP VHEAACRQF RDLPPVARKDRQFWSLWWVDYFCHSHFLWHVVFVVLGVVGHYRAVLQMSRIV WLDAGRAF</p>
Source	in vitro E.coli expression system
Target Names	IZH3
Protein Names	Recommended name: ADIPOR-like receptor IZH3
Expression Region	1-419
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.