



Recombinant Human Probable G-protein coupled receptor 63 (GPR63)

Product Code	CSB-CF871626HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q9BZJ6
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MVFSAVLTAFHTGTSNTTFVVYENTYMNITLPPPFQHPDLSPLLRYSFETMAPT GLSSLT VNSTAVPTTPAAFKSLNLPLQITLSAIMIFILFVSFLGNLVVCLMVYQKAAMRSAI NILL ASLAFADMLLAVLNMPFALVTILTTRWIFGKFFCRVSAMFFWLFVIEGVAILLIISI DRF LIIVQRQDKLNPYRAKVLIASWATSFCVAFPLAVGNPDLQIPSRAPQCVFGYT TNPGYQ AYVILISLISFFIPFLVILYSFMGILNTRLRHNALRIHSYPEGICLSQASKLGLMSLQR PF QMSIDMGFKTRAFTTILILFAVFIVCWAPFTTYSLVATFSKHFYQHNFFEISTW LLWLC YLKSALNPLIYYWRIKKFHDACLDMMPKSFKFLPQLPGHTKRRIRPSAVYVCGE HRTVV
Source	in vitro E.coli expression system
Target Names	GPR63
Protein Names	Recommended name: Probable G-protein coupled receptor 63 Alternative name(s): PSP24-2 PSP24-beta
Expression Region	1-419
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.