



# Recombinant Pan troglodytes Melanin-concentrating hormone receptor 1 (MCHR1)

<b>Product Code</b>	CSB-CF013584EQV
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	Q5IJ49
<b>Product Type</b>	Transmembrane Protein
<b>Immunogen Species</b>	Pan troglodytes (Chimpanzee)
<b>Sequence</b>	MSVGAMKKGVGRAVGLGGGSGCQATEEDPLPNCGACAPGQGGRRWRLPQ PAWVEGSSARL WEQATGTGWMDLEASLLPTGPNASNTSDGPDNLTSAGSPPRRTGSISYINIIMP SVFGTIC LLGIIGNSTVIFAVVKKSKLHWCNNVPDIFIINLSVVDLLFLLGMPFMIHQLMGNG VWHF GETMCTLITAMDANSQFTSTYILTAMAIDRYLATVHPISSTKFRKPSVATLVICLL WALS FISITPVWLYARLIPFPGGAVGCGIRLPNPDLDLYWFTLYQFFLAFALPFVVITAA YVRI LQRMTSSVAPASQRSIRLRTKRVRTAIAICLVFFVCWAPYYVLQLTQLSISRPT LTFVY LYNAAISLGYANSCLNPFVYIVLCETFRKRLVLSVKPAAQGQLRAVSNAQTADE ERTESK GT
<b>Source</b>	in vitro E.coli expression system
<b>Target Names</b>	MCHR1
<b>Protein Names</b>	Recommended name: Melanin-concentrating hormone receptor 1 Short name= MCH receptor 1 Short name= MCH-R1 Short name= MCHR-1 Alternative name(s): G-protein coupled receptor 24 MCH-1R Short name= MCH1R Short name=
<b>Expression Region</b>	1-422
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.