



Recombinant Rabbit Secretin receptor (SCTR)

Product Code	CSB-CF020865RB
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O46502
Product Type	Transmembrane Protein
Immunogen Species	Oryctolagus cuniculus (Rabbit)
Sequence	HLVGAPPRLCDVLWVLQEERDQCLQEELERERLGEEQPVPGCQGLWDNVSC WPSSAPGRMV ELECPRFLRMLTNSNGSLFRNCTQDGWTETFPRLACGVSMNDSSHERQH AYLLKLVKVM YTVGYSSSLVMLLVALGILCAFRRHLHCTRNYIHMHLFLSFILRALSNIKDAVLFS SDDA IHCDHRVGGCKLVMVFFQYICIMANYAWLLVEGLYLHSLLVVSFFSERKCLQGF VVLGWGS PAMFVTSWAVTRHFLEDSGCWDINANAIIWWVIRGPVILSILINIFILFINILRILTR KLR TQETRGQDMNHYKRLARSTLLLIPLFGVHYIVFVFSPEGAMEIQLFFELALGSF QGLVVA VLYCFLNGEVQLEVQKKWQQWHLWEPPLCPVALSSSFNNGTSSLNSTKACPS GRSRDTCK VSII
Source	in vitro E.coli expression system
Target Names	SCTR
Protein Names	Recommended name: Secretin receptor Short name= SCT-R
Expression Region	22-445
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.