



Recombinant Human Neuromedin-U receptor 1 (NMUR1)

| | |
|--------------------------|---|
| Product Code | CSB-CF884509HU |
| Storage | Store at -20°C, for extended storage, conserve at -20°C or -80°C. |
| Uniprot No. | Q9HB89 |
| Product Type | Transmembrane Protein |
| Immunogen Species | Homo sapiens (Human) |
| Sequence | <p>MTPLCLNCSVLPGLDLYPGGARNPMACNGSAARGHFDPEDLNLTDEALRLKYL GPQQTELF MPICATYLLIFVVGAVGNGLTCLVILRHKAMRTPTNYLFLSLAVSDLLVLLVGLP LELYE MWHNYPFLLVGGCYFRLLFEMVCLASVLNVTALSVERYVAVVHPLQARSM VTRAHVRR VLGAVWGLAMLCSLPNTSLHGIRQLHVPCRGVPDSAVCMLVRPRALYNMVV QTTALLFF CLPMAIMSVLYLLIGLRLRRERLLLMEAKGRGSAARSRYTCRLQQHDRGRR QVTKMLF VLVVVFVFGICWAPFHADRVMWSVVSQWTDGLHLAFQHVHVISGIFFYLGSAANP VLYSLMS SRFRETFQEALCLGACCHRLRPRHSSHLSRMTTGSTLCDVGSVLSGWSVHPLA GNDGPEAQ QETDPS</p> |
| Source | in vitro E.coli expression system |
| Target Names | NMUR1 |
| Protein Names | Recommended name: Neuromedin-U receptor 1 Short name= NMU-R1 Alternative name(s): G-protein coupled receptor 66 G-protein coupled receptor FM-3 |
| Expression Region | 1-426 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | full length protein |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |