



# Recombinant Pig Vasoactive intestinal polypeptide receptor 1 (VIPR1)

<b>Product Code</b>	CSB-CF637928PI
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	Q28992
<b>Product Type</b>	Transmembrane Protein
<b>Immunogen Species</b>	Sus scrofa (Pig)
<b>Sequence</b>	<p>GLQQEEDY LQMIKVQHKQC LEEAQL ENETSGCSKMWDNLTCWPATPRGQV  VVLACPLIF  KLFSP TQGLN VSRNCTDEGWTPLEPGPYIACGMDDKASGLDEQQTVFYNSV  KTGYTIGY  SLSLAALLVATAILSLFRKLHCTRNYIHMHLFISFILRATAVFIKDLALFDSEESDH  CSK  GSVGCKAAV VLFQYCVMANFFWLLVEGLYLHTLLAVSFFSERKYFWGYIFVG  WGV PSTFI  M VWTVVRIHFEDYGCWDTIHSSLWWI KAPILASILVNFILFIRIIGILVQKLRPPD  VGK  SDNSPYSRLAKSTLLLIPLFGVHYIMFAFFPDNFKA EVKMFELIVGSFQGCVVA  ILYCF  LNGEVQAELRRKWR RWHQQGV LGWDSKYQHPSGG SNGDTCSTQV SMLTR  VSPSARRSSSF QAEVSLV</p>
<b>Source</b>	in vitro E.coli expression system
<b>Target Names</b>	VIPR1
<b>Protein Names</b>	Recommended name: Vasoactive intestinal polypeptide receptor 1 Short name=VIP-R-1 Alternative name(s): Pituitary adenylate cyclase-activating polypeptide type II receptor Short name= PACAP type II receptor Short name= PACAP-R
<b>Expression Region</b>	32-458
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.