



Recombinant Rabbit Cholecystokinin receptor type A (CCKAR)

Product Code	CSB-CF004772RB
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O97772
Product Type	Transmembrane Protein
Immunogen Species	Oryctolagus cuniculus (Rabbit)
Sequence	MDAVASLLGNASGIPPPCELGLDNETLFCLDQPPPSKEWQPAVQILLYSLIFLL SVLGNT LVITVLIRNKRMRVTNIFLLSLAISDLMLCLFCMPFNLIPNLLKDFIFGSALCKTT TYL MGTSVSVSTLNLVAISLERYGAICKPLQSRVWQTKSHALKVIAATWCLSFAIMT PYPIYS NLVPFTKTNNQTANMCRFLLPSDVMQQAWhTFLLILFLIPGIVMMVAYGMISL ELYQGI KFDASQKKSARKERKASTGSGRFEDNDGCYLQRSKPTRQLELQQLSGGGGGR VSRIRSSSS AATLMAKKRVIRMLMVIVVLFFLCWMPIFSANAWRAYDTVSAERRLSGTPISFIL LLSYT SSCVNPIIYCFMNRFRRLGFMATFPCCPNPGPPGPRAEAGEEEEEGRTRASLS RYSYSHM SASAPPS
Source	in vitro E.coli expression system
Target Names	CCKAR
Protein Names	Recommended name: Cholecystokinin receptor type A Short name= CCK-A receptor Short name= CCK-AR Alternative name(s): Cholecystokinin-1 receptor Short name= CCK1-R
Expression Region	1-427
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.