



# Recombinant Dog Cholecystokinin receptor type A (CCKAR)

<b>Product Code</b>	CSB-CF688568DO
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	Q5D0K2
<b>Product Type</b>	Transmembrane Protein
<b>Immunogen Species</b>	Canis lupus familiaris (Dog) (Canis familiaris)
<b>Sequence</b>	MEVADSLGNGSDVPPPCELGLENETLVCLEQPRAAKEWQPAVQILLYSLIFLL SVLGNT LVITVLIRNKRMRVTNIFLLSLAVSDLMLCLFCMPFNLIPNLLKDFIFGSAVCKTT TYF MGTSVSVSTFNLVAISLERYGAICKPLQSRVWQTKSHALKVIATTWCLSFTIMT PYPIYS NLVPFTKTNNQTANMCRFLLPNDVMQQSWHTFLLLILFLIPGIVMMVAYGLISLE LYQGI KFDAIQKKSARDRNPSTGSSGRYEDGDGCYLQKARPRRRLRQLSTPGSGR LNRIRSTS STANLMAKKRVIRMLMIVVLFFLCWMPIFSANAWRAYDTASAERRLSGTPISF ILLLSY TSSCVNPIIYCFMKNRFRFLGFLATFPCCPHPGPPGPRGEVGEEEEGRTTGASL SRYSYSH MSASAPGP
<b>Source</b>	in vitro E.coli expression system
<b>Target Names</b>	CCKAR
<b>Protein Names</b>	Recommended name: Cholecystokinin receptor type A Short name= CCK-A receptor Short name= CCK-AR Alternative name(s): Cholecystokinin-1 receptor Short name= CCK1-R
<b>Expression Region</b>	1-428
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.