



# Recombinant Meleagris gallopavo Beta-4C adrenergic receptor (ADRB4C)

<b>Product Code</b>	CSB-CF337330MQV
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P43141
<b>Product Type</b>	Transmembrane Protein
<b>Immunogen Species</b>	Meleagris gallopavo (Wild turkey)
<b>Sequence</b>	<p>MTPLPAGNGSVPNCSSWAAVLSRQWAVGAALSITILVIVAGNLLVIVIAIAKTPRLQ  TMTNV  FVTSLACADLVMGLLVVPPGATILLSGHWPYGTVVCELWTSLDVLCVTASIEL  CAIAVD  RYLAITAPLQYEALVTKGRAWAVVCMVWAI SAFISFLPIMNHWWRDGADEQAV  RCYDDPR  CCDFVTNMTY AIVSSTVSFYVPLLVMIFVYVRVFAVATRHVQLIGKDKVRFLQE  NPSLSS  RGGRWRRPSRLLAIKEHKALKTLGIIMGTFTLCWLPFFVANIIVFCRPLVPDQL  FLFLN  WLGYNVNSAFNPIIYCRSPDFRSFRKLLCCPRRADRRLLHAAPQDPQHCSAFC  PRGDPME  DSKAVDPGHLREDSEVQSGRREENASSHGGGHQQRPLGECWLQGMQSM  LCEQLDEFTST EMPAGPSV</p>
<b>Source</b>	in vitro E.coli expression system
<b>Target Names</b>	ADRB4C
<b>Protein Names</b>	Recommended name: Beta-4C adrenergic receptor Alternative name(s): Beta-4C adrenoreceptor Short name= Beta-4C adrenoceptor
<b>Expression Region</b>	1-428
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.