



Recombinant Dictyostelium discoideum Frizzled/smoothened-like sans CRD protein G (fscG)

Product Code	CSB-CF687544DKK
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q54DR6
Product Type	Transmembrane Protein
Immunogen Species	Dictyostelium discoideum (Slime mold)
Sequence	QSLPSLPSPITYKSSQCRGDLVYRNVSNRMEDEKIGFTFVQYNGTYQSCVIPC PSPFFTL DEWNKFLYMSLVMGTISFLCGLFLLITYSPIVNKTHNRHTIGVMCMSFGVCLAM CSDMWN FGSNFTDQKSICPSPGQYLTTNSRCLGSGIVLQFGGVFGFLNWTLLSFDLFM NIKGIIT KNYDKYYFVATFIIAIIFTFVPIVNDQYSMSYIGLGCWLGSVAVYQLIFFWILLSICLI VS SVFIILILKEIYIIIKQSKQKTSKGNIRPLLCITVTSFAFFYMFFYYISIVIEGDYER ILNEYTDCLMDPTKDVSECKFPRMSVANEFVLLCLRLLGIGAFIFYGINKEVKKI WLNS FWFNNSFVGKYIGSKRSMGNDITNSYASKAYSKNYNNNNSINSYNSGLELSIID MSCNKD DNFKPIIK
Source	in vitro E.coli expression system
Target Names	fscG
Protein Names	Recommended name: Frizzled/smoothened-like sans CRD protein G
Expression Region	25-453
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.