



Recombinant Escherichia coli Uncharacterized protein ynbD (ynbD)

Product Code	CSB-CF301179ENV
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P76093
Product Type	Transmembrane Protein
Immunogen Species	Escherichia coli (strain K12)
Sequence	<p>MLQGAGWLLLLAPFFFFTYGSLNQFTAVQDLNSHDIPSQVFGWETAIPFLPWTI VPYWSL DLYGFSLFVCSTTFEQRRLLVHRLILATVMACCGFLLYPLKFSFIRPEVSGVTG WLFSQL ELFDLPYNQSPSLHIILCWLLWRHFRQHLAERWRKVCGGWFLIAISTLTTWQ HHFIDVI TGLAVGMLIDWMVPVDRRWNYQKPDQRRRIKIALPYVVGAGSCIVLMELMMMI QLWWSVWL CWPVLSLLIIGRGYGGLGAIITGKDSQGKLPVAVYWLTLPCRIGMWLSMRWFC RRLEPVS KMTAGVYLGAFPRHIPAQNAVLDTFVPRGRATKDRLYFCVPMMLDLVVPEEG ELRQAVA MLETLREEQGSVLVHCALGLSRSALVVAWLLCYGHCKTVNEAISYIRARRPQI VLTDEH KAMLRLWENR</p>
Source	in vitro E.coli expression system
Target Names	ynbD
Protein Names	Recommended name: Uncharacterized protein ynbD
Expression Region	1-430
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.