



Recombinant Meleagris gallopavo D (2) dopamine receptor

Product Code	CSB-CF529149MQV
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O73810
Product Type	Transmembrane Protein
Immunogen Species	Meleagris gallopavo (Wild turkey)
Sequence	MDPLNLSWYNTGDRNWSEPVNESSADQKPQYNYAVLLTLLIFVIVFGNVLVC MAVSREK ALQTTTNYLIVSLAVADLLVATLVMPWVVYLEVVGWRFHSRIHCDIFVTLDMMM CTASIL NLCAISIDRYTAAAMPMLYNTRYSSKRRVTVMACVWVLSFAISSPILFGLNKAD ERECI IANPAFVVYSSVVSFYVPFIVTLLVYVQIYMVLRRRRKRHTKRSSHGLSDTHA PLKDKC THPENVKLGTIVVKSNGSFQVNRKCEAESHKMEMEMMSSTSPPERTIVKAA APSNHQL VVPVASSRSTLDSPGKVEKNGHAKENLHTAKVFEIQSMPNGKTRSTLLKAMNR RKLSQQK EKKATQMLAIVLGVFIICWLPFFITHILNMHCDCNIPPAMYSFTWLGYN SAVN PIIYT TFNIEFRKAFMKILHC
Source	in vitro E.coli expression system
Target Names	DRD2
Protein Names	Recommended name: D(2) dopamine receptor Alternative name(s): Dopamine D2 receptor
Expression Region	1-436
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.