



Recombinant Meleagris gallopavo Vasoactive intestinal polypeptide receptor (VIPR1)

Product Code	CSB-CF841943MQV
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q91085
Product Type	Transmembrane Protein
Immunogen Species	Meleagris gallopavo (Wild turkey)
Sequence	GLVVEVEVWWWRRWRFGGGGCIMLEIEEERSQCPAEITEDNQTSGCRRQWD NITCWPEAQV GAVVVKPCPKYFRLLTTF LGNVS RNCTS QGWTDVYPAPYAVACGYDSTAHQG KEQTAFYG TVKTGYTIGHTLSLIALTAAMIILCLFRKLHCTRNYIHMHLFMSFIMRAIAVFIKDV TLF ESGEPEHCFVSSVGCKAMMVFFQYCVMANFFWLLVEGLYLHTLLVISFFSERK YFWWYIL IGWGAPSVFITAWTVVRIYFFNVGCWEEIIESPIWWIKTPILVSILVNFILFICIIRIL VQKLHSPDVGHNETSQYSRLAKSTLLLIPLFGIHYIMFAFFPDNFKAQVKLVFEL VVGFSF QGFVVAVLYCFLNGEVQAELKRKWRRWHLERFLGSDMKYHHP SLG SNGTNF STQISMLTK CSPKTRRCSSFQAEFSLV
Source	in vitro E.coli expression system
Target Names	VIPR1
Protein Names	Recommended name: Vasoactive intestinal polypeptide receptor Short name= VIP receptor Short name= VIP-R
Expression Region	20-457
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.