



Recombinant Rat Histamine H3 receptor (Hrh3)

Product Code	CSB-CF886187RA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q9QYN8
Product Type	Transmembrane Protein
Immunogen Species	Rattus norvegicus (Rat)
Sequence	MERAPPDGLMNASGTLAGEAAAAGGARGFSAAWTAVLAALMALLIVATVLGN ALVMLAFV ADSSLRTQNNFFLLNLAISDFLVGAFCIPLYVPYVLTGRWTFGRGLCKLWLVVD YLLCAS SVFNIVLISYDRFLSVTRAVSYRAQQGDTRRAVRKMALVWVLAFLLYGPAILSW EYLSGG SSIEGHCYAEFFYNWYFLITASTLEFFTPFLSVTFFNLSIYLNQRRLRLDGG REAG PEPPDAQPSPPPAPPSCWGCWPKGHGEAMPLHRYGVGEAGPGVEAGEAA LGGGSGGGAA ASPTSSSGSSSRGTERPRSLKRGSKPSASSASLEKRMKMVSQSITQRFRLSR DKKVAKSL AIIVSIFGLCWAPYTLLMIIRAACHGRCIPDYWYETSFWLLWANSVNPVLYPLC HYSFR RAFTKLLCPQKLKVQPHGSLEQCWK
Source	in vitro E.coli expression system
Target Names	Hrh3
Protein Names	Recommended name: Histamine H3 receptor Short name= H3R Short name= HH3R
Expression Region	1-445
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.