



Recombinant Danio rerio Calcitonin gene-related peptide type 1 receptor (calcrla)

Product Code	CSB-CF731565DIL
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q68EK2
Product Type	Transmembrane Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Sequence	<p>SPEVNESQQQHPQNVYHDIGVTRNKIVTAQFECYQKIMKDNSQDRRGPVCNR TWDGWLCW DDTEAGITSEQHCPDYFQDFDPTMVTKICTESGWFLHPESNRTWTNFTRC NLHTTEGR RTAMNLFYLALIGHGLSLTSLFISLGIFFHFKSLSCQRITLHKNLFFSFVLNSIITII WL TAVANNQELVQQNPISCKISQFIHLYIFGCNYFWMLCEGIYLHTLIVVAVFAEKQ HLMWY YLLGWGFPLIPATIHAVARSYYYNDNCWISSNTSLLYIIHGPICAAMLVNLFFLLN IVRV LITKLVTHQAKSSLYMKAVRATLILVPLLGIQYVLLPYKPSGRVSAEIYDYIMHIL MHY QGLLVATIFCFFNGEVQAVLRRHWNQYRIQFGSTITQSDALRSASYTASSITEV QGCYSI DGHTEHLNGKNYHDFDNAIKPENPFA</p>
Source	in vitro E.coli expression system
Target Names	calcrla
Protein Names	Recommended name: Calcitonin gene-related peptide type 1 receptor Short name= CGRP type 1 receptor Alternative name(s): Calcitonin receptor-like receptor
Expression Region	24-470
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.