



Recombinant Mouse Probable G-protein coupled receptor 61 (Gpr61)

Product Code	CSB-CF803973MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q8C010
Product Type	Transmembrane Protein
Immunogen Species	Mus musculus (Mouse)
Sequence	MESSPIQSSGNSSTLGRALQTPGPSTASGVPELGLRDVASESVALFFMLLLD LTAVAGN AAVMAVIAKTPALRKVFVVFHLCLVDLLAALTMPLAMLSSSALFDHALFGEVA CRLYLF LSVCFVSLAILSVSAINVERYYYVHPMRYEVRMTLGLVASVLVGVVVKALAM ASVPVLG RVYWEEGAPSVNPGCSLQWSHSAYCQLFVVVFAVLYFLLPLILIFVYCSMFR VARVAAM QHGPLPTWMETPRQRSESLSSRSTMVTSSGAHQTTPHRTFGGGKAAVLLA VGGQFLLCW LPYFSFHLYVALSAQPISAGQVENVVTWIGYFCFTSNPFFYGCLNRQIRGELSK QFVCFE KAAPEEELRLPSREGSIEENFLQFLQGTSENWVSRPLPSPKREPPPVD FRIP GQIAEET SEFLEQQLTSDIIMSDSYLRPAPSPRLES
Source	in vitro E.coli expression system
Target Names	Gpr61
Protein Names	Recommended name: Probable G-protein coupled receptor 61
Expression Region	1-449
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.